

# The Animas® Vibe® Insulin Pump and CGM\* System

## Factory Defaults

Setup Screen	Factory Defaults	Options
<b>Time</b>		
Hour	Mode: 12 Hr	12 or 24 Hr
<b>Sounds</b>		
A - Audio Bolus	H (High)	Vib (Vibrate), L (Low), M (Medium), H
N - Normal Bolus	H	OFF, Vib, L, M, H
T - Temp Basal	OFF	OFF, Vib, L, M, H
Alert	H	OFF, Vib, L, M, H
Reminder	H	Vib, L, M, H
Warning	H	Vib, L, M, H
Alarms	H	Vib, L, M, H
<b>Advanced</b>		
<b>Setup ADV 1 (Ability to program 12 segments)</b>		
I:C ratio	12:00 AM, 1:15	1:2 – 1:160
ISF	12:00 AM, 1:50	1:5, 1:600
BG Target	12:00 AM, 120 +/- 10	60 – 250 mg/dL, (+/-) 0-50 mg/dL
<b>Setup ADV 2 BOLUS</b>		
Adv. Bolus	OFF	OFF, ON
Reminders	OFF	OFF, ON
Delivery	NRML	SLOW, NRML
Basal Programs	1	1 or 4
<b>Setup ADV 3 Max LIMITS</b>		
Max Basal	5.00 U/Hr	0.00 U/Hr – 25.00 U/Hr
Max Bolus	16.00 U	0.00 U – 35.00 U
Max Daily	100 U	1 U – 600 U
2 hours	50 U	1 – 100 U
<b>Setup ADV 4</b>		
Language	English	English, alternate local language
Display Timeout	30 secs	15, 30, 45, or 60 seconds
Contrast	8	1 – 10
Battery	Lith	Lithium or alkaline
<b>Setup ADV 5 AUTO OFF</b>		
Enable	OFF	OFF, ON
Suspend Duration	12 Hr	1 Hr – 24 Hr
<b>Setup ADV 6</b>		
Low Cartridge Warning	20 units	10, 20, 30, 40, 50
Occlusion Sensitivity	L	Low or High
<b>Setup ADV 7 Audio Bolus</b>		
Enable	OFF	OFF, ON
Delivery Step	1.0 U	0.1 U, 0.5 U, 1.0 U, 5.0 U
<b>Setup ADV 8 INSULIN on BOARD</b>		
IOB 2	OFF	OFF, ON
Duration	4.0 Hr	1.5 Hr – 6.5 Hr

\*Continuous glucose monitoring.

# The Animas® Vibe® Insulin Pump and CGM System

## Factory Defaults (cont.)

Setup Screen	Factory Defaults	Options
<b>Setup ADV 9 Sick Days</b>		
BG over ____mg/dL	240	60 – 400 mg/dL
Check Ketones every ____Hrs	4 Hrs	1 Hr – 24 Hrs
Check BG every ____Hrs	2 Hrs	1 Hr – 10 Hrs
<b>CGM Warning Sounds</b>		
High Alert	H	Vib, L, M, H
Low Alert	H	Vib, L, M, H
Rise Rate	H	Vib, L, M, H
Fall Rate	H	Vib, L, M, H
Range	H	Vib, L, M, H
Other	H	Vib, L, M, H
<b>CGM Setup</b>		
High Alert	200 mg/dL	120 – 400 mg/dL
High Alert Snooze	0 (OFF)	0 – 300 minutes
High Alert Enable	YES	YES, NO
Low Alert	80 mg/dL	60 – 100 mg/dl
Low Alert Snooze	0 (OFF)	0 – 300 minutes
Low Alert Enable	YES	YES, NO
Rise Rate	3 mg/dL/min	2, 3 mg/dL/min
Rise Rate Enable	YES	YES, NO
Fall Rate	3 mg/dL/min	2, 3 mg/dL/min
Fall Rate Enable	YES	YES, NO
Out of Range Snooze	30 mins	21 mins – 201 mins
Out of Range Enable	YES	YES, NO
Transmitter ID	00000	60000 – 7YYYY

### Important Safety Information

The Animas® Vibe® Insulin Pump and CGM System is intended for the delivery of insulin and for continuous glucose monitoring (CGM) for the management of insulin-requiring diabetes. The Animas® Vibe® System's CGM, which includes the Dexcom G4® PLATINUM Sensor and Transmitter, is indicated for detecting trends and tracking patterns in persons age 2 and older. The system is intended for single patient use and requires a prescription.

**Contraindications:** Insulin pump therapy is not recommended for people unwilling or unable to test their blood glucose four to six times per day, unwilling or unable to see their healthcare professional regularly, or whose vision or hearing does not allow recognition of pump alerts, warnings, and alarms. The Animas® Vibe® Insulin Pump must be removed before MRI or CT scan, and the Dexcom G4® PLATINUM Sensor and Transmitter must be removed before MRI, CT scan, or diathermy treatment. Taking acetaminophen-containing medications while wearing the sensor may falsely raise sensor glucose readings.

**Warnings:** Connecting to the pump before receiving the necessary training could result in serious injury or death. Younger children may inadvertently press the pump buttons and deliver insulin, which can lead to hypoglycemic events. Caregivers are responsible for helping to ensure safe and effective delivery of insulin to people in their care, including using the safety features on the Animas® Vibe® System to help prevent injury. CGM glucose readings are not to be used to make treatment decisions. The sensor and transmitter do not replace a blood glucose meter. The sensor and transmitter are not to be used during pregnancy or while on dialysis. Sensor placement is only approved for sites under the skin of the belly (abdomen) in adults and the belly or upper buttocks for ages 2 to 17 (pediatrics). Patients should seek professional medical help if a sensor breaks and no portion of it is visible above the skin, or if there is infection or inflammation. Any broken sensors or adverse events should be reported to Customer Service. Contact Customer Service at 1-877-937-7867 or visit [www.animas.com/importantsafetyinformation](http://www.animas.com/importantsafetyinformation) for detailed indications for use and safety information.