



Transferring settings to your new Animas® insulin pump

Congratulations on receiving your new Animas® insulin pump. This document is designed to assist you with transferring your current insulin pump settings to your new Animas® insulin pump. Record your personalized settings from your current Animas® insulin pump in the spaces provided. Once recorded, you must manually enter these settings into your new Animas® insulin pump.

Contact your healthcare team if you have any questions about the settings you are entering. Failure to use settings determined by your healthcare team will likely result in incorrect insulin delivery amounts.

YOUR CURRENT ANIMAS® INSULIN PUMP SETTINGS

Your current basal program rates are located in your Animas® insulin pump in [Main Menu > Basal](#). An "A" will appear to the left of the basal program that is currently active. Enter your current basal settings below.

Basal Program Rates							
Program 1		Program 2*		Program 3*		Program 4*	
Weekday		Other		Weekend		Exercise	
Total U		Total U		Total U		Total U	
Start	U/HR	Start	U/HR	Start	U/HR	Start	U/HR
12:00 A		12:00 A		12:00 A		12:00 A	

*You may have 4 basal programs displayed in the Basal menu if you activated multiple basal programs in Set-up > Advanced 2 (Basal Programs).

If you require assistance obtaining settings from your insulin pump or programming these settings, please call the phone number that is located on the back of the insulin pump.

Enter your setting below or check where appropriate.

Main Menu > Bolus > Reminders†

Reminder 1 ON OFF
Time AM PM

Reminder 2 ON OFF
Time AM PM

BG Check ON OFF
After Bolus 1 Hr 2 Hr 3 Hr 4 Hr

Main Menu > Setup > Time/Date

Ensure that the time, date and year that you set in your insulin pump is accurate. If you select the 12-hour time format (Mode), be sure the desired AM or PM selection is correctly displayed.

Main Menu > Setup > Sound

High (H), Medium (M), Low (L), Vibrate(V), or OFF

N-Bolus (Normal Bolus)
A-Bolus (Audio Bolus)
T-Basal (Temporary Basal)
Alert
Reminder
Warning
Alarm

Main Menu > Setup > Advanced 1 Audio Bolus

Enable ON OFF
Delivery Step 0.1U 0.5U 1.0U 5.0U

Main Menu > Setup > Advanced 2

Bolus

Adv. Bolus ON OFF
Reminders ON OFF
Delivery NRML SLOW

Basal

Programs 1 4

Main Menu > Setup > Advanced 3 Max Limits

Max Basal U / Hr
Max Bolus U
Max Daily U



Animas Corporation

965 Chesterbrook Blvd

Wayne, PA 19087

Phone: 1-877-YES-PUMP (1-877-937-7867)

Email: comments@animas.com

www.animas.com

†You may have your personal Reminders feature turned on or off. The Reminders feature is located in the Main Menu > Setup > Advanced 2 (Reminders).

When this feature is turned on, your insulin pump will display Reminders in the Bolus Menu.

‡ Advanced 7, 8, 9a, and 9b in 1200+ and 1250 models. In 1200 model, ISF and BG Target settings are in the ezBolus menu and no range (+/-) setting is available.

§ Advanced 10 in 1200+ and 1250 models.

|| Advanced 11 in 1200+ and 1250 models. Sick Day screen not available on 1200 model.

All product names and trademarks are the property of their respective owners.

© 2015 Animas Corporation 2015/06 42047400 Rev A, ANM 15-4452A Rx Only.