

My Insulin Pump Settings



Basal Programs and Basal Rates: Main Menu > Basal

Basal Program 1 Weekday		Basal Program 2 Other		Basal Program 3 Weekend		Basal Program 4 Exercise	
Basal Programs 2-4 are available if you have opted for 4 Basal Programs (Main Menu > Setup)							
Start Time	Units/Hr	Start Time	Units/Hr	Start Time	Units/Hr	Start Time	Units/Hr
12:00A		12:00A		12:00A		12:00A	

Personal Settings: Main Menu > Setup

Insulin to Carb (I:C) Ratio		Insulin Sensitivity Factor (ISF)		BG Target/Range		
Start Time	I:C (gm)	Start Time	ISF (mg/dL)	Start Time	mg/dL	+/- mg/dL
12:00A	1 Unit:	12:00A	1 Unit:	12:00A		

Other Pump Settings: Main Menu > Setup

Sound – High, Med, Low, Vib, OFF

N-Bolus (Normal)	A-Bolus (Audio)	T-Basal (Temp Basal)	Alert	Reminder	Warning	Alarm

Bolus Options

Adv bolus On Off		Contrast	Reminders On Off		Delivery Normal Slow		Basal Programs 1 or 4	
Max Limits	Display Time Out		Auto-Off Alarm	Low Cart Warning	Occlusion Sensitivity	Insulin On Board (IOB)	Audio Bolus	
Basal	15 sec	(1 – 10)	Enable: OFF ON When Enabled suspend in ___ Hrs	10 units 20 units 30 units 40 units 50 units	L (Low) H (High)	ON OFF Duration ___ Hrs	Enable: OFF ON Delivery Step: ___ U	
Bolus	30 sec							
Daily	45 sec							
2Hr	60sec							

Reminder Settings: Main Menu > Bolus

Reminder 1			Reminder 2			Reminder 3		
OFF	ON	Time = _____	OFF	ON	Time = _____	OFF	ON	After bolus _____ Hr