



To Whom It May Concern:

_____ has diabetes and wears an Animas® Insulin Pump to control it. As the manufacturer of this pump, Animas recommends that the pump **does not** go through:

- X-ray screening that is used for carry-on or checked luggage
- Whole Body Imaging Technology

X-ray technologies may damage the pump, resulting in a potential impact on insulin delivery to the patient and subsequent adverse event. We recommend that the pump be screened using alternate methods.

Additionally, this patient:

- will be wearing an infusion set with a small catheter placed under the skin
- may be wearing a Dexcom® Continuous Glucose Monitor Sensor and Transmitter also placed under the skin.

Neither of these two items should be removed from the body. These items can be subjected to a metal detector.

This patient may have diabetes supplies packed in carry-on luggage for immediate access and safety, including:

- a back-up insulin pump, should the first pump fail while traveling
- insulin pump cartridges
- infusion sets
- insulin
- insulin syringes
- blood glucose monitoring supplies: meter, test strips, lancing device, lancets
- quick acting glucose tablets or gel

If you have any questions regarding the information above, please contact Animas Corporation at 877-767-7373 (United States and Canada) or +001-610-644-8990 (outside the United States and Canada).

Sincerely,

Animas Corporation
965 Chesterbrook Boulevard
Wayne, PA